

primera[™]

TENS/NMES WITH HAN WAVEFORM



Multi-functional TENS and NMES Device

- Extensive range of programs – 8 TENS and 6 NMES Programs
- The HAN TENS Waveform: A unique dispense and disperse mode of stimulation which modulates between low (2 Hz) and high (70 Hz) frequencies.
- TENS programs used for temporary pain relief
- NMES programs stimulate healthy muscles in order to facilitate and improve muscle performance.
- Complete package includes: PRIMERA Unit, Carrying Case, Battery, Leadwires, Electrodes, Electrode Positioning Back Belt and User's Manual
- Two-year Warranty



primera[™]

TENS/NMES WITH HAN WAVEFORM

Technical Specifications

- **Dual channel:** individually isolated circuits
- **Type:** Constant Current
- **Low Battery Indicator:** If the battery goes below 6.9 volts +/- 0.2 Volts the battery symbol will flash on/off once every second
- **Waveform:** Asymmetrical, rectangular bi-phasic with zero DC current
- **Open Electrode Detect:** If an open circuit is detected at the output of channel A or B the output current will be reset at zero
- **TENS and NMES Amplitude:** 0 - 80 mA
- **Programs Include:** Continuous TENS, Burst TENS, Modulated TENS, HAN TENS, and NMES

ORDERING INFORMATION

77630* PRIMERA[™] TENS/NMES Kit (Over-the Counter Device)
77615 PRIMERA[™] TENS/NMES Kit (Prescription Device)

42181 Dura-Stick Plus 2" x 3.5" Rectangle Electrodes (Includes 10 packs of 4 electrodes in each pack)
42183 Dura-Stick Plus 2" Square Electrodes (Includes 10 packs of 4 electrodes in each pack)

* Primera (PN 77630) is an "Over-the-Counter (OTC)" device and is not reimbursed by private or government payors including Medicare

Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this brochure do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.



DJO, LLC | A DJO Global Company

T 800.592.7329 D 800.242.8329

1430 Decision Street | Vista, CA 92081-8553 | U.S.A.

DJOglobal.com/chattanooga