



## Transitional 172 Cervical Collar

### Optimal bruk

Halskrave som begrenser fleksjon, ekstensjon og rotasjon av halsryggraden. ProCare® Transitional 172 halskrave ble utformet for å gi optimal støtte og komfort. Kraven begrenser bevegelsemulighetene kraftig, samtidig som den gir reduserte trykkpunkter som kan føre til hudmasesasjon og/eller sår. De myke, pustende skumputene gir økt komfort og oppmuntrer til økt pasientcompliance. Avtakbare, åpne skumputeceller fjerner fuktighet fra huden for god hygiene og komfort samt minimalt med hudirritasjoner. Røntgen-, CT- og MRI-gjennomskinnelig.

### Fordeler og egenskaper

- Effektiv begrensning av bevegelse
- Redusert belastning av huden
- Støtte som kan tilpasses
- Utbyttbare puter
- Utbyttbare polstringer

### Bestillingsinformasjon

TRANSITIONAL 172 CERVICAL COLLAR (ADULT)		
Artikkelnr.	Beskrivelse	Størrelse
79-83283	Transitional 172 Collar, Adult Short	Short
79-83285	Transitional 172 Collar, Adult Regular	Regular
79-83287	Transitional 172 Collar, Adult Tall	Tall
79-83288	Transitional 172 Collar, Adult X-Tall	X-Tall
79-83363	Transitional 172 Collar Set*, Adult Short	Short
79-83365	Transitional 172 Collar Set*, Adult Regular	Regular
79-83367	Transitional 172 Collar Set*, Adult Tall	Tall
79-83368	Transitional 172 Collar Set*, Adult X-Tall	X-Tall

TRANSITIONAL 172 CERVICAL COLLAR (PEDIATRIC)		
Artikkelnr.	Beskrivelse	Størrelse
79-83272	Transitional 172 Collar, Ped-1	Ped-1
79-83273	Transitional 172 Collar, Ped-2	Ped-2
79-83275	Transitional 172 Collar, Ped-3	Ped-3
79-83277	Transitional 172 Collar, Ped-4	Ped-4
79-83278	Transitional 172 Collar, Ped-5	Ped-5
79-83351	Transitional 172 Collar Set*, Ped-1	Ped-1
79-83352	Transitional 172 Collar Set*, Ped-2	Ped-2
79-83353	Transitional 172 Collar Set*, Ped-3	Ped-3
79-83354	Transitional 172 Collar Set*, Ped-4	Ped-4
79-83355	Transitional 172 Collar Set*, Ped-5	Ped-5

